



# BREAKTHROUGH ADVANCES IN SYNCHROTRON XRPD:

a Powerful Product-Characterization Technique with  
Extended Applications to Pharmaceuticals

Presented by

**Arturo Araque, PhD**  
Head of Excelsus Consortium

**Fabia Gozzo, PhD**  
Head of Excelsus Structural Solutions

Produced by:

 **webinar**

Excelsus



# Excelsus

SCIENTIFIC ENGINEERING

## Quality Assurance

## Validation

Commissioning

Qualification




Process Validation

Cleaning Validation









## Where in pharmaceutical production?

- **Drug characterization**
- **Product purity**
- **Raw materials control**
- **Counterfeiting prevention**
- **API fingerprinting**
- **Cross-contaminant quantification**

# From Research tool to Commercial use

- 1) Validate**
- 2) Ensure Quality**
- 3) Communicate**



From Research tool to Commercial use – Validation

**Commissioning**

**Installation Qualification**

**Operational Qualification**

**Performance Qualification**

**Methods Validation**

From Research tool to Commercial use – Quality

**Standard Operating Procedures**


**Bill of Materials**

**Certificate of Analysis**

**License documentation**

# From Research tool to Commercial use - Announce

- 1. The advantages of using synchrotron over conventional laboratory XRPD**
- 2. Technical instrumental advances**
- 3. Benefits of this new technology**



## Outlook

- Role of structural analysis in the pharmaceutical industry
- X-Ray Powder Diffraction and advantages of synchrotron radiation
- Synchrotron XRPD applications: indexing & structural solutions, radiation damage control, kinetic studies, Quantitative Phase Analysis (L.o.D, L.o.Q)
- Easy and affordable access to state-of-the-art synchrotron XRPD through Excelsus

# Why is structural analysis of drugs relevant to the pharmaceutical industry?


- Polymorphism of pharmaceuticals and their **structure ↔ properties**
- Microstructural properties (e.g. influence of stress and strain, particle size and domain)




# Polymorphism of pharmaceuticals

Polymorphism is the ability of substances to crystallize in solid state phases according to different arrangements or conformations of the molecules in the crystal lattice


Example of  
Bupivacaine Hydrochloride



**Form B at 112 °C**




**Form D at 20 °C**



Gozzo, Masciocchi, Griesser, Niederwanger, 2010. Personal Communication

# Why is polymorphism such an important property?



**Characterization and control of pharmaceutical properties through characterization and control of their solid-state structure**

## Properties that are influenced by the solid state structure of the substance and, therefore, influenced by polymorphism:

- Solubility
- Pharmacokinetics and pharmacodynamics
- Thermodynamic properties (e.g. stability of drugs) → **in-situ non-ambient studies**
- Spectroscopic properties
- Mechanical properties (e.g. hardness, compressibility, tableting, tensile strength)

### ICH Guideline (Q6A):

*Differences in these [polymorphic] forms could, in some cases, **affect the quality or performance of the new drug products**. In cases where differences exist which have been shown to affect drug product performance, bioavailability or stability, then the appropriate solid state should be specified (3.3.1. New Drug Substance).*

*Identification testing should establish the identity of the new drug substance(s) in the new drug product and should be able to **discriminate between compounds of closely related structure** which are likely to be present (3.2.2. New Drug Product)*

***New analytical technologies**, and modifications to existing technology, are continually being developed. Such technologies should be used when justified (ICH scope)*

## The famous case of Ritonavir\*

- Abbott's protease inhibitor Norvir, generically Ritonavir (AIDS treatment)
- Not bioavailable from the solid state → Marketed in 1996 as a oral-liquid and semisolid capsules (both formulations contained ethanol/water based)
- **No crystal form control** performed before distribution on the market because the ICH guidelines stated that *“For a drug product that is a solution, there is little scientific rationale for polymorph control”*
- In Summer 1998, Norvir semisolid capsules supplies were threatened as a result of its spontaneous irrepressible interconversion into a new thermodynamically stable form precipitating out of the semisolid formulated product, due to unexpected supersaturation and nucleation conditions.
- The new form had much lower solubility with greatly reduced bioavailability
- Drug **out of market for 1 year** until Ritonavir polymorphism was fully understood,
- The metastable form I could then be stabilized with a new gel capsule formulation.

 **During development, characterization of all possible crystal forms is needed, even in liquid formulations**

*\*Ritonavir: an Extraordinary Example of Conformational Polymorphism*  
Bauer et al, Pharmaceutical Research 18(6). 2001, 859-866

# Polymorphic studies play a key role throughout the whole life-cycle of products

## Compound selection

- Identification and characterization of individual polymorphic forms
- Selection of desired form based on therapeutic properties



## Technical development

- Development of DS and DP manufacturing processes to ensure high and reproducible content of desired polymorphic form via fingerprinting
- Polymorphic studies for impurity detection and stability studies
- Consider crystal engineering (e.g. co-crystallization\*) to enhance yield and stability of desired polymorphic form



## Commercial production

- Polymorphic characterization to support (1) process validation, (2) comparability studies following process changes, and (3) investigations to assess impact of deviation on product quality

Synchrotron  
Radiation X-ray  
Powder Diffraction  
is a unique and  
powerful technique  
for such studies

\* Example of carbamazepine: rational co-crystallization with saccharin that prevents the formation of hydrated forms that have strongly reduced solubility!  
(see *Organic Crystal Engineering*, Eds. Tiekink, Vittal & Zaworotko, Wiley 2010)



## Additional areas where polymorphic studies play a key role

- ★ Fight against **counterfeit drugs**:
  - ➡ Trademark infringements
  - ➡ Serious attempts to human health → no or wrong API content, hazardous adulterants, substituted and not labeled ingredients


*See: “The fight against counterfeit drugs”* (<http://www.cbsnews.com/video/watch/?id=7359537n>)

- ★ **Intellectual Property**: → scientific support to litigations involving patent issues

# **Synchrotron X-Ray Powder Diffraction**

## **(SR-XRPD)**

# X-Ray diffraction by crystals





Position of the diffracted peaks → size and dimension of the unit cell

Intensity ratios of the diffracted peaks → type and location of atoms in the unit cell


Full Width at Half Maximum (FWHM) →  
of the diffracted peaks

Intrinsic properties of the materials  
(i.e. microstructural analysis)

# X-Ray diffraction from a polycrystalline or a powder sample



$$|\mathbf{r}^*| = \frac{1}{d_{hkl}}$$




Debye-Scherrer (1916)

In a powder diffraction pattern the 3D of the reciprocal space collapse onto 1D



the information in the orientation space gets lost



Extracting detailed information from a powder pattern requires very high-quality data  
→ the **collimation** and **monochromaticity** of a **synchrotron X-ray beam** allow for a great improvement of the **angular resolution** of the acquired patterns, whereas the **high brilliance** of a synchrotron source drastically reduces measurements times allowing the study of the **kinetic of transformations** in fraction of seconds




# Why powder diffraction?

- Single crystal not always available
- The majority of APIs and excipients in pharmaceutical substances and products are powders (polycrystalline or amorphous powders)
- Quantitative Phase Analysis of mixtures (e.g. drug product)
- In-situ, non-ambient XRD: powders often a better/the only choice
- Radiation damage → XRPD parallel 1D detection
- Total scattering techniques


# Synchrotron X-Ray Powder Diffraction

**What makes it a powerful analytical tool?**

## Synchrotron Radiation and beamline optics




## State-of-the-art diffractometers




Swiss Light Source-Materials Science beamline  
Powder Diffraction station

## Outstanding detection systems



Multicrystal Analyser

Schmitt et al, 2003,  
Bergamaschi, Schmitt et al, 2010



MYTHEN II


## SR Properties:

- High Spectral Brightness:  
[photons/s/0.1%  
bandwidth/mm<sup>2</sup>/mrad<sup>2</sup>] → 10<sup>12</sup>-10<sup>15</sup>  
ph/sec in small beams (μm to mm)
- Tunable photon energy
- Polarization
- Time structure
- Coherence

Resolution: 1 arcsec  
Accuracy: ±2 arcsec  
Precision: ±1 arcsec

**Highest flexibility to accommodate  
all kinds of sample environments**

## Synchrotron Radiation and beamline optics




## State-of-the-art diffractometers



Swiss Light Source-Materials Science beamline  
Powder Diffraction station

## Outstanding detection systems



Multicrystal Analyser

## SR Properties:

- High Spectral Brightness:  
[photons/s/0.1%  
bandwidth/mm<sup>2</sup>/mrad<sup>2</sup>] → 10<sup>12</sup>-10<sup>15</sup>  
ph/sec in small beams (μm to mm)
- Tunable photon energy
- Polarization
- Time structure
- Coherence


Resolution: 1 arcsec  
Accuracy: ±2 arcsec  
Precision: ±1 arcsec

**Highest flexibility to accommodate  
all kinds of sample environments**

- ❖ Ultra-high resolution (better than 0.003°)  
independent of sample dimension
- ❖ Angular selection of diffracted beam and  
all related advantages:
- ❖ Fluorescence suppression
- ❖ Independence of resolution on sample size  
and sample positioning
- ❖ Independence on transparency effect
- ❖ High S/N and S/B
- ❖ Long measurements (15 min-several hours)



## Synchrotron Radiation and beamline optics



## State-of-the-art diffractometers



Swiss Light Source-Materials Science beamline  
Powder Diffraction station

## Outstanding detection systems


- ❖ Solid state microstrip modular detector
  - ❖ 1D position sensitive
  - ❖ 120° angular coverage
  - ❖ Single photon counting read out → fluorescence-free bkg
  - ❖ Large dynamic range (24 bits)
  - ❖ 0.004° inherent angular resolution
    - FWHM ~0.009° for 0.2 mm capillaries
    - ~0.022° for 0.5 mm capillaries
    - ~0.044° for 1.0 mm capillaries
  - ❖ High d-spacing resolution ( $d_{\min}=0.24 \text{ \AA}$  at 28keV)
- Time resolved powder diffraction:**
- ❖ Average acquisition time 0.5-1s, down to 100msec or nsec (pump-and-probe)
  - ❖ The acquisition can run 30000 faster than using the single channel crystal analyzer

## SR Properties:

- High Spectral Brightness:  
[photons/s/0.1%  
bandwidth/mm<sup>2</sup>/mrad<sup>2</sup>] → 10<sup>12</sup>-10<sup>15</sup>  
ph/sec in small beams (μm to mm)
- Tunable photon energy
- Polarization
- Time structure
- Coherence

Resolution: 1 arcsec  
Accuracy: ±2 arcsec  
Precision: ±1 arcsec


**Highest flexibility to accommodate  
all kinds of sample environments**



MYTHEN II



## Bupivacaine Hydrochloride - form D



- 1 mm capillary,
- Mythen data at 50% reduced intensity
- No radiation damage up to 3min




Large counting statistics  
in subsec acquisition times →

In-situ kinetic studies of organic  
compounds!

Gozzo F. , 2008


### Experimental details:

- MYTHEN I and II
- In-situ measurements in the STOE furnace
- Several cycles to check the reversibility of form B



**DSC measurements**  
(Courtesy U. Griesser,  
University Innsbruck)

### In-situ T-dependent powder diffraction patterns




Niederwanger, V., Gozzo, F. & Griesser, U. (2009). J. Pharm. Sci. 98, 1064–1074

Gozzo F., Masciocchi N., Griesser U. & Niederwanger, V. 2010, private communication


### Experimental details:

- MYTHEN I and II
- In-situ measurements in the STOE furnace
- Several cycles to check the reversibility of form B



**DSC measurements**  
(Courtesy U. Griesser,  
University Innsbruck)

### In-situ T-dependent powder diffraction patterns




Niederwanger, V., Gozzo, F. & Griesser, U. (2009). J. Pharm. Sci. 98, 1064–1074

Gozzo F., Masciocchi N., Griesser U. & Niederwanger, V. 2010, private communication


### Experimental details:

- MYTHEN I and II
- In-situ measurements in the STOE furnace
- Several cycles to check the reversibility of form B



**DSC measurements**  
(Courtesy U. Griesser,  
University Innsbruck)

### In-situ T-dependent powder diffraction patterns




Niederwanger, V., Gozzo, F. & Griesser, U. (2009). J. Pharm. Sci. 98, 1064–1074

Gozzo F., Masciocchi N., Griesser U. & Niederwanger, V. 2010, private communication


### Experimental details:

- MYTHEN I and II
- In-situ measurements in the STOE furnace
- Several cycles to check the reversibility of form B



**DSC measurements**  
(Courtesy U. Griesser,  
University Innsbruck)

### In-situ T-dependent powder diffraction patterns




Niederwanger, V., Gozzo, F. & Griesser, U. (2009). J. Pharm. Sci. 98, 1064–1074

Gozzo F., Masciocchi N., Griesser U. & Niederwanger, V. 2010, private communication

# Structural solution with Simulating Annealing and Direct Methods of the RT phase of the Bupivacaine Hydrochloride (Form A)


2 seconds XRPD acquisition time



Details of sa: 2 fragments, Cl free,  
 7 torsions, 15 DoF

# Structural solution with **Simulating Annealing** and **Direct Methods** of the RT phase of the Bupivacaine Hydrochloride (Form A)

2 seconds XRPD acquisition time



Details of sa: 2 fragments, Cl free,  
 7 torsions, 15 DoF



Example of structural solution of a subtle case of polymorphism:

## CARPROFEN


Non-steroidal anti-inflammatory drug (NSAIDs)

Identification of two forms:

- Re-crystallization from hot solution of commercial sample in ethanol (*Form I*)
- High-Temperature *Form II* as solid-solid transition

Thermal, spectroscopic, lab-XRPD and *ab-initio* SR-XRPD and post-experimental DFT calculations

Low-T measurements and high-pressure investigations



Outcome:

Form I → Form II subtle **isostructural transformation**


Form I and II are **conformational polymorphs**

Evidence for **configurational disorder**

**SR-XRPD successful in identifying two very close polymorphs that could not be distinguished with laboratory XRPD data**

G. Bruni, F. Gozzo et al. J. Pharm. Sci. 2011, **100**(6), 2321-2332

10 sec acquisition time with optimized optics,  
reduced intensity beam and MYTHEN II



### Form I

Monoclinic  $P2_1/c$

$a=8.2855(3)$

$b=5.7071(3)$

$c=27.087(1)$

$\beta=96.901(5)$

$Z=4$

$R_{wp}(\%)=1.77$ ;  $GoF=2.85$

### Form II

Monoclinic  $P2_1/c$

$a=8.5563(3)$


$b=5.7923(3)$

$c=26.778(1)$


$\beta=100.473(5)$

$Z=4$

$R_{wp}(\%)=2.09$ ;  $GoF=3.26$



Form I




Form II




- *In-situ* hydrogen absorption at 15 bar
- *In-situ* desorption by connecting the cell to a vacuum pump
- Continuous measurements using the  $\mu$ strip detector while the reaction takes place
- Acquisition times between 5 and 20 sec per pattern, depending on the reaction kinetic.

Courtesy of Prof. Cerny, Uni. Geneva - Joubert et al. Acta Materialia, **54** (2006), 713-719


### hydrogen uptake



### diffraction pattern



### phase content



New beamline optics+ Mythen II → 5-10 faster

t = 40 s

# Quantitative Phase Analysis (QPA)

QPA refers to the ability of quantitatively state the abundance of the different phases that constitute a mixture.

## Why is this relevant?

- ❑ **Polymorphic purity**: detect and quantify unwanted polymorphic forms in both drug substance and drug product
  - Level of Detection (LoD)
  - Level of Quantitation (LoQ)
- ❑ Assess the **polymorphic composition** in drug substance and product
- ❑ In formulated materials, the **API/excipients relative proportion** is paramount and needs to be kept under control
- ❑ **Degree of Crystallinity** in amorphous/crystalline mixtures

## QPA analytical methods

Several are the analytical methods used to obtain quantitative phase related information:

- Spectroscopic methods (e.g. Raman and Infrared spectroscopy, Mass spectroscopy, Nuclear Magnetic Resonance spectroscopy)
- Thermal Methods (e.g. Differential Scanning Calorimetry, ThermoGravimetric Analysis)
- Diffraction Methods → XRPD Direct method

**More on QPA:** *Quantitative Phase Analysis*, **Madsen & Scarlett**, Chapter II in “Powder Diffraction: Theory and Practice”, **Dinnebier & Billinge** Editors, 2008, The Royal Society of Chemistry, Cambridge, UK

*Introduction to X-Ray Powder Diffractometry*, **Jenkins & Snyder**, Publisher, Wiley 1996

*Modern Powder Diffraction*, **Bish & Post Editors**, Reviews in Mineralogy, Vol.20, 1989.


*Quantitative phase analysis by combining the Rietveld and the whole-pattern decomposition methods*, **Giannini, Guagliardi & Millini**, *J.Appl.Cryst* (2002). 35, 481.

## Single-peak methods:


- Intensity ratio  $\frac{I_{unknown}}{I_{standard}}$  of one or several diffraction peaks
- No need of structural information but prone to systematic errors caused by preferential orientation (PO) and reflection overlapping

## Whole-pattern methods:

- Use of full diffraction patterns
- minimization of systematic errors, peak overlap not a problem
- Preferential Orientation (PO) can be modeled
- Accuracy close to X-Ray fluorescence elemental analysis, with the advantage of being sensitive to structural differences → direct QPA of polymorphs

- 
- Rietveld Method (Rietveld, JAC 1969. 2, 65) → QPA from the structural model
  - Partial Or No Known Crystal Structure (PONKCS) Method (Madsen & Scarlett, 2006)
    - needed pure phases (or almost pure) and *ad-hoc* mixtures pure phases+standard in known % wt; PONKCS works well for quantitation of amorphous phase!
  - QUANTO/Wilson-plot technique (Giannini, Guagliardi et al, 2002)
    - needed pure phases, cell parameters and SG, but NO NEED of *ad-hoc* mixtures

# Standard 99.9% Haloperidol + 0.1% Indomethacin



All synthetic physical mixtures provided by G. Bruni, Uni. Pavia, Italy

The same mixture at 3 different locations on the capillary:

1. Halo/Indom:  
99.82%/0.18%
2. Halo/Indom:  
99.91%/0.09%
3. Halo/Indom:  
99.97%/0.03%

Average: Halo/Indom:  
99.9%/0.1%


Refinement on merged data:

Halo/Indom:  
99.91%/0.09%




Laboratory XRPD  
1.0 mm capillary; transmission geometry  
Fast Inel detector (1024 sec acquisition time)

Lab data  
Courtesy Ian Madsen



Lab and SR XRPD data collected on identical synthetic mixtures



- Polymorphic forms may have a significant impact on the **quality** or **performance** of pharmaceutical products
- A careful characterization of the polymorphism of drug substance and drug product should therefore play a key role throughout the **whole life-cycle** of pharmaceutical products
- Polymorphic studies have also started to play a key role during **patent litigations** and in the fight against **counterfeit drugs**
- **Synchrotron-Radiation Powder Diffraction** has become a unique and very powerful tool for polymorphic studies, such as kinetic analyses, the identification of closely related polymorphic forms and high-sensitivity quantitative phase analyses
- This use is in line with the regulatory expectations (ICH guidelines and FDA guidance) that newly available analytical technologies are used for **continuous improvements** in process understanding and product characterization

### Level I:

You internally analyze your synchrotron data.

- we assist you in the choice of the appropriate experimental conditions;
- we prepare the samples
- organize and perform the synchrotron measurements
- deliver the data to you with a written report


### Level II:

In addition to Level I services

- we analyze and interpret your SR-XRPD data
- provide a written report

### Level III:

- In addition to Level I and II services,
- we provide scientific consulting
- work in close collaboration with you
- assist you in the design of new experiments



**Company web site:**  
**[www.excels.us](http://www.excels.us)**

Contact us via email to discuss your specific needs:  
[fabia.gozzo@excels.us](mailto:fabia.gozzo@excels.us) or +32 470 990851  
[arturo.araque@excels.us](mailto:arturo.araque@excels.us) or +1 510 5937540

➤ **Paul Scherrer Institute Direction** and its **Technology Transfer Department**

➤ **Swiss Light Source** and **Materials Science** beamline:

Ch. Quitmann, F. van der Veen

Phil Willmott, Antonio Cervellino, Nicola Casati, Michael Lange, Dominik Meister

Bernd Schmitt, Anna Bergamaschi (PSI detector group)

➤ Gianna Bruni (mixtures for QPA scientific project), University of Pavia, Italy

➤ Ian Madsen and Nikki Scarlett (QPA scientific project), CSIRO, Victoria, AU

➤ Cinzia Giannini & Brunella Aresta (QPA scientific project), CNR, Bari, Italy

➤ Roberto Bellotti (QPA scientific project), University of Bari, Italy

➤ Irene Margiolaki (powder diffraction of proteins), University of Patras, Greece

➤ Scientists and colleagues always open to scientific discussions :

Radovan Cerny, Carmelo Giacobazzo, Angela Altomare, Anna Moliterni, Rosanna Rizzi,

Michele Saviano, Thomas Laube, Arnt Kern, Bernd Hinrichsen, Peter Stephens, Thomas

Hartmann, Piero Macchi, Norberto Masciocchi, Antonella Guagliardi, Meitian Wang,

Martin Fuchs, ...